

HOT PLATE



// CABINET FEATURES:

- 1020 carbon steel;
- Corrosion-resistant;
- Epoxy powder textured coated steel;
- Rubber feet;
- Safety fuse holder;
- Stainless steel support between cabinet and plate;
- Polycarbonate frontal label.

// INCLUDED ACCESSORIES:

- Safety fuse.

// TECHNICAL DATA:

TEMPERATURE RANGE	From 50°C up to 300°C (other temperatures, under analysis)
PLATE MATERIAL	Brushed aluminium with thermal insulation
TEMPERATURE CONTROL	Electronic adjustment 50°C in 50°C
PRE SET ON/OFF	Luminous Indication
POWER CABLE	3 pins - Brazilian standard (NBR14136)
HEATING SYSTEM	Shielded resistance - easy to replace
POWER SUPPLY	110V or 220V - to be defined

// MODEL DETAILS:

MODEL	TOP PLATE SIZE	OVERALL DIMENSIONS	TOTAL WEIGHT (KG)	POWER
SSCA	25x30 cm	15x25x30 cm	10	600 W
SSCA	30x40 cm	15x30x40 cm	12	1100 W
SSCA	35x45 cm	15x35x45 cm	15	2000 W



// OPTIONAL:

- Stainless steel hot plate.

Notice: Any product presented can be produced according to customer's needs. All models come with English language instruction manual. 1-year warranty against manufacturing defects.