

AIR STERILIZER



// CABINET FEATURES:

- 1020 Carbon Steel;
- Anti-corrosion treatment;
- Electrostatic textured epoxy coating;
- Wall Mounting fittings;
- Safety fuse;
- Easy-to-install.

// INCLUDED ACCESSORY:


- Power cable 3-pin plug (NBR14136);
- 10 filters;
- Wall mounting support;
- Screws and bushings;
- Safety fuse.

// TECHNICAL DATA:

AIR CIRCULATION	Forced ventilation, low noise
LAMPS	Germicidal UVC Fluorescent - 18W
FILTER	Particulate filter for easy cleaning
POWER CABLE	3-pin plug (NBR 14136)
POWER	110v or 220V

// MODEL DETAILS:

MODEL	CAPACITY	DIMENSIONS (HEIGHTxWidthxDepth)	LAMPS	POWER
STERIAIR 120	120 m ³ /h	84x25x8 cm	4 UNITS	60 W



// OPTIONAL:

- Outside in 430 stainless steel;
- Hour meter.

Notice: Any product presented can be produced according to customer's needs. All models come with English language instruction manual. 1-year warranty against manufacturing defects.